

# Development of a Scintillator/Tungsten Electromagnetic Calorimeter with Silicon Photomultipliers

Joseph Proulx

Elliot Smith

Keith Drake

University of Colorado at Boulder


# The Colorado Group

Shirley Choi, Keith Drake, Joshua Elliot, Eric Erdos, Chris Geraci, Jack Gill, Jason Gray, Kyle Miller, Sarah Moll, Uriel Nauenberg, Gleb Oleinik, Joseph Proulx, Elliot Smith, Jiaxin Yu, Jinlong Zhang

# Outline

- Our design
- WLS fiber test
- Silicon photomultiplier studies
- Detector module studies with cosmic rays
- Future work and upgrades

# Our design - basic overview




- 1/2-radiation-length Tungsten plate absorber (1.75mm)
- 2mm-thick scintillator sensitive material
- 5cm square scintillator tiles, offset by 2.5cm in each layer
- WLS fiber light guides
- SiPM detectors – easy to calibrate, small, low-voltage/low power, inexpensive


# Our design - possible variations

- Tyvek vs. radiant mirror
- Orientation of WLS fibers
- Mounting SiPM's on tiles
- Ganging multiple tiles together


# Testing our design



- WLS fiber stability




- SiPM performance



- sensitivity of tiles

# WLS fiber long-term test



We have been observing the possible degradation of the WLS fiber over time for the last year. We have studied fibers bent at a radius of 2cm and 6cm, annealed and not annealed.



# WLS fiber long-term test

## Annealing process

- heat molds to 104C
- slide fibers into mold
- Fibers cool into shape








# WLS fiber long-term test

## Results

- One year of data shows no significant degradation of fibers
- Uncertainty dominated by variation in LED output




# SiPM Studies




● Single-photon sensitivity

# SiPM Studies




● Single-photon sensitivity




● Temperature Dependence


# SiPM Studies



● Single-photon sensitivity



● Temperature Dependence





● Time resolution

# SiPM Studies


## Basic readout scheme

- Amplify signal from SiPM by 200
- Separate signal into 4 channels, each with a 2.5ns delay
- Read each channel at 10ns intervals
- Merge data from 4 channels




# SiPM Studies - setup


- We drive a single SiPM with an ultraviolet LED coupled to WLS fiber with neutral-density filters to control light input



# SiPM Studies - setup




# SiPM Studies - initial results



- Early results demonstrated ability to sense very low light levels, but did not show resolution on the single-photon level





# SiPM Studies - cold setup




# SiPM Studies - low-temperature results

- Cold-temperature tests demonstrated improved signal for similar noise (but still no single-photon resolution)





# SiPM Studies - low-temperature results




# SiPM Studies

- New devices have much more stable response and better gain




# SiPM Studies - Single-photon separation!




# SiPM Studies - time resolution

- Reading pulse in 10ns slices, we can still separate photons




# SiPM Studies - time resolution

- Reading pulse in 10ns slices, we can still separate photons



# SiPM Studies - time resolution


- Reading pulse in 10ns slices, we can still separate photons






# SiPM Studies - time resolution

- Reading pulse in 10ns slices, we can still separate photons



# SiPM Studies - time resolution

- Reading pulse in 10ns slices, we can still separate photons



# SiPM Studies - time resolution


- Photons from LED arrive over a 40ns time interval
- Response of SiPM is about 10ns
- Performance is similar in each of the SiPM's we studied

# Detector module tests


We use cosmic rays to test our calorimeter module.

- Particles ionize scintillator in triggers and module tile
- Triggers deposit light in PMT's
- Coincidence in PMT signals used as trigger for data collection
- Signal from tile converted by SiPM; data collection same as for SiPM tests

# Detector module tests




# Detector module tests




# Detector module tests

We observe a signal of several photons over 100ns



# Detector module tests


Integration of the entire pulse results in noise that fades the single-photon resolution






# Detector module tests

- Reading pulse in 10ns slices, we can again separate photons




# Detector module tests

- Reading pulse in 10ns slices, we can again separate photons




# Detector module tests

- Reading pulse in 10ns slices, we can again separate photons




# Detector module tests

- Reading pulse in 10ns slices, we can again separate photons




# Detector module tests

- Reading pulse in 10ns slices, we can again separate photons




# Detector module tests

- Reading pulse in 10ns slices, we can again separate photons



# Detector module tests

- Reading pulse in 10ns slices, we can again separate photons



# Future Plans

- Improve data analysis
  - Better peak finding
  - Account for pedestal voltage
  - In principle, this can be solved by better electronics
- Study uniformity of tile response
  - Need code to calibrate SiPM response
- Investigate feasibility of ganging tiles
  - Light guides?
  - Larger SiPM surface area?
- Incorporate findings into ongoing simulation studies